

The Wymo Wheelchair Hoist

The Wymo Wheelchair Hoist provides independence to a motorist who requires a wheelchair in their daily life.



The Wymo Wheelchair Hoist is an Australian product

How it works?

The hook of the hoist is lowered & hooked onto the wheelchair by the driver. The driver uses the switch box to gently guide the wheelchair up onto the hoist from the frame on the car roof.

The driver is now free to drive anywhere, anytime. Once the driver reaches their destination, the Down switch is pressed on the control box, which smoothly delivers the wheelchair to the ground.

- ✓ Operation can be from either the Driver or Passenger side of the vehicle.
- ✓ Can be fitted to most roof rack systems.
- ✓ Light weight (17kg)
- ✓ Powered by the car's own battery – either through the cigar lighter socket or direct wired.



Frank's Engineering Pty Ltd

135 Bakers Road North Coburg Vic 3058

Ph: 03 9354 0400 | Fax: 03 9354 0344

Email for Quotes: sales@franksengineering.com.au

www.franksengineering.com.au